Your
**Port Access
Advantage***

BardPort*
Implanted Ports
Portales Implantados

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Table Of Contents

INTRODUCTION

Introduction	1
What is a Bard implanted port?	2
How is the port placed?	2
How the port is used	3

ACCESS

How your port is accessed	4
-------------------------------------	---

CARE

Care of your port	4
-----------------------------	---

HEPARIN LOCKS

About heparin locks	5
-------------------------------	---

POWER INJECTION

Feel the New Standard of Care*	6
Identification	6
Patient Discharge Packet	7
Tell your clinician	7
CECT scans / power injection	8

Power

Tell Your Clinician.	9
------------------------------	---

Q & A

Common questions and answers	10
--	----

Introduction

Your doctor has recommended an implanted port from Bard Access Systems for your IV therapy because your treatments require frequent delivery of medications directly into your bloodstream. With your implanted port, you'll be able to receive IV therapy without having to undergo repeated needlesticks in your peripheral (arm or leg) veins.


Please read all of the information contained in this Patient Guide. It is designed to give you the knowledge to understand and feel comfortable with your implanted port.

This information is provided to you as a reference. Your doctor and nurse are still your best sources of information. Any questions or concerns that you have should be directed to them. Always follow their directions concerning the use and maintenance of your implanted port. You can also visit www.portadvantage.com for more information about implanted ports from Bard Access Systems.

What is a Bard implanted port?

Your Bard implanted port is a small device placed completely beneath your skin in a short procedure. It is a cylinder with a hollow space inside that is sealed by a soft top called a *septum*, which is made of a self-sealing rubber material. The port device connects to a small, flexible tube called a catheter that is inserted inside one of the large central veins that deliver blood to your heart. When a special needle is put into the soft top of the device, it creates “access” to your bloodstream, meaning that medications and fluids can be given and blood samples withdrawn.

Implanted ports are available in a variety of shapes and sizes. Your physician has chosen the Bard implanted port which is most appropriate for your body size and the type of therapy you are receiving.

How is the port placed?

The Bard implanted port is placed during a brief procedure that usually is preformed under local anesthesia. You may even be allowed to return home the same day. Your doctor will place the port just beneath the skin and insert the end of the catheter into the blood vessel that has been selected to receive the medication. Your doctor will locate the port in a position most suitable for your treatment. Frequently, the port is placed on the upper chest in an area just below the collarbone. The port will be visible merely as a small “raised area” beneath the skin.

How the port is used


Your implanted port allows clinicians to easily deliver medications or fluids or withdraw blood samples without having to repeatedly stick your peripheral veins directly with a needle. This makes it more comfortable for you. Because the port device places medications into the large central veins instead of the small peripheral veins, the medications mix more thoroughly in the blood, diluting them so they are less harmful to your vascular system.

The port can be used to deliver medication in two ways. One method, called “bolus injection”, delivers medication all at once with the needle left in place only for a short time. The other method is a slow delivery of medication called a “continuous infusion”. For this type of medication, a dressing will usually be placed over the needle to hold it in the port for a longer period of time and a smaller tube will connect it to the medication source.

The venous port also provides an easy method for collecting blood samples. When a blood sample is needed, a special needle is placed through the soft top of the port device which connects to the catheter.

How your port is accessed

Your clinicians will use the implanted port whenever they need to infuse medications or fluids into your body or withdraw samples of your blood. To do this, they will first access the port by inserting a special needle into the soft top of the device. For most patients, there is only a mild pricking sensation felt during needle insertion. Frequently, the sensation of the needle insertion decreases over time. Ask the clinician about what to expect in your own procedure.

Care of your port

During the first couple of days after receiving your port, it will be important to avoid heavy exertion and to follow any special instructions given to you by your doctor for care of the small incision. Once the incision has healed, you will not have to take any special care of the port and you can resume your normal daily activities.

About heparin locks

Sometimes when blood clots, it can block the catheter, preventing medications and fluids from flowing through it. However, blood will not clot when it is thinned with a medication called “heparin.” To help prevent clots from forming, implanted ports with open-ended catheters are typically filled with sterile heparinized saline after each use. This process is called a “heparin lock.” If your implanted port will not be used for long periods of time, the clinician will typically change the heparin lock every four weeks.



If you have been diagnosed as being allergic to heparin or as having Heparin Induced Thrombocytopenia (HIT), it is important to remind the clinician of this fact anytime your device is accessed.

Feel the New Standard of Care*

The PowerPort* Implantable Port is a new kind of vascular access device for your IV therapy treatments. It offers the unique ability to provide access for power-injected Contrast-Enhanced Computed Tomography (CECT) scans. Power-injected CECT scans produce superior images of your body to help the medical team better manage your care. With a PowerPort* implanted port, you'll be able to receive IV therapy and CECT scans without having to undergo repeated needlesticks in your arm or leg veins.

Identification

Several innovative features help to identify the PowerPort* device. Along with its distinctive purple color and unique triangular shape, the PowerPort* implanted port features a proprietary triangular arrangement of three bumps on the top of the septum called Palpation Points.

Patient Discharge Packet

To further enhance recognition, patients also receive a special Patient Discharge Packet that includes an Identification Card, Bracelet and Key Ring when they are implanted with the device.

Tell your clinician


If you have a PowerPort* implanted port, a trained medical professional can identify the unique shape and design that distinguishes it from traditional ports. Tell your clinician that you have been implanted with the PowerPort* device beforehand whenever your port is accessed for a procedure.

CECT scans / power injection

Contrast-Enhanced Computed Tomography (CECT) scans are simple, safe and non-invasive procedures that provide quick and accurate information for disease diagnosis and treatment. These scans are many times more sensitive and radiologists can distinguish small differences in your soft tissues that may not be detected with conventional x-rays.

Before performing a CECT scan, the radiology team will inject a special fluid in your body to help produce clearer pictures during the CECT scan procedure. For best results, the fluid, called contrast media, is infused at a high rate into your bloodstream. This process is called power injection.

A PowerPort* Implantable Port used with a PowerLoc* Safety Infusion Set has the unique ability to allow clinicians to perform power-injected CECT scans without having to make a needlestick in your arm or leg veins.

Tell Your Clinician

As a new patient with an Implanted Port Device, the best way you can be involved is to share information and concerns with your clinicians. Speak up! Ask questions about anything that concerns you or if you notice anything that seems unusual.

What to Report to Your Clinician

Here are some important things to share with your clinician:

- You have an Implanted Port.
- If you have a **Bard PowerPort*** device.
- If you notice any redness or inflammation at the site of your implanted port after the incision heals
- If you have a fever.
- If you have allergies to any medications or materials.
- If you have an allergy to heparin.
- If you have heparin induced thrombocytopenia (HIT)
- If you have ever been prescribed anticoagulant (blood thinning) medications such as heparin or warfarin.
- If you have been previously treated with radiation.
- If you have ever been diagnosed with lung disease.
- If you have ever been diagnosed with, or treated for, venous thrombosis.
- If you have ever been diagnosed with any tissue disease or suffered from tissue erosion.
- If you have ever been diagnosed with or tested for “pinch off” syndrome.
- If other clinicians have ever had difficulty withdrawing blood or infusing fluids through your implanted port, including whether other clinicians have ever required you to change position to allow blood or fluid to flow.

Common questions and answers

Q: Will the implanted port affect my daily activities?

A: Once the incision heals following implantation, you should be able to return to your normal daily activities, such as bathing, swimming or jogging. Ask your doctor or nurse about specific activities and the appropriate time to resume them.

Q: Will I need to wear a bandage over the implanted port?

A: A bandage will be required until your incision heals. After the incision has healed, a bandage is not required when the port device is not being used. If you are receiving continuous infusion of fluids, a bandage may be applied to stabilize and protect the needle while it is in place.

Q: Do I have to stop wearing certain types of clothing?

A: Ask your doctor or nurse because the answer will depend on where your port device is placed.

Q: Who pays for the implanted port?

A: Insurance policies vary, so check with your insurance company first.

Q: Will the implanted port activate security alarms?

A: Security systems may detect the small amount of metal in the device. If it does occur, simply show your device Patient Identification Card.

Q: How long will I have my implanted port?

A: The implanted port can stay in place for as long as your doctor determines that you need it.

Q: Will my port need to be accessed when not in use?

A: Yes, it will need to be flushed every 4 weeks.

Q: Can the device be removed if I no longer need it?

A: Yes. When no longer needed, the implanted port can be removed in a simple procedure similar to the one used to implant it.

Q: Can I get a CT procedure with an implanted PowerPort* device?

A: Yes, The materials used in the PowerPort* system are non-ferromagnetic and are safe for use in CT and CECT procedures.

Q: Can I get an MRI procedure with an implanted PowerPort* device?

- A:** The device can be scanned safely under the following conditions: Non-clinical testing has demonstrated the device is MR Conditional. It can be scanned safely under:
- static magnetic field of 3 Tesla or less
 - spatial gradient field of 330 Gauss/cm or less
 - maximum specific absorption rate (SAR) of 6 W/kg for 30 minutes of scanning.

In non-clinical testing, the device produced a temperature rise of less than 0.5 °C at a maximum specific absorption rate (SAR) of 6 W/kg for 30 minutes of MR scanning in a 3T Siemens Trio with software version VA25.

For Minimal Image artifact

MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the device. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic implant.

Q: What do I do if I have a question or suspect a problem?

- A:** Remember that following the procedure, you will experience some redness and tenderness in the area of the incision which normally subsides in the first 24 to 48 hours. However, contact your doctor or nurse if you notice unusual changes in the skin area over the port such as increased swelling, redness or soreness. If you experience pain, fever, chills, shortness of breath or dizziness, contact your doctor immediately.

Q: What do I do if my doctor or nurse is unavailable?

- A:** You will be provided with a wallet-sized Medical Alert Identification Card that contains important information about the port. With this information, other qualified medical personnel will be able to provide assistance.

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Warning: An issued or revision date for these instructions is included for the user's information. In the event two years have elapsed between this date and product use, the user should contact Bard Access Systems to see if additional product information is available.

Revised Date: August, 2008

*Bard, BardPort, PowerPort, PowerLoc, X-Port, "Feel the New Standard of Care", "Your Port Access Advantage" and the color purple are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate. ChronoFlex is a registered trademark of CardioTech International, Inc.

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MC-0115-00
0808R

Acknowledgement Card

I, _____, have received
print name of patient or care giver
the booklet "**Bard Access Systems** Implanted Port Patient
Information" from _____.

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(Firma del paciente)

(Fecha)

You should carry your Medical Alert
Identification Card with you at all times.

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Caution: Use Only Non-Coring
Needles

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