



## Device Identification Web Page Synchromed<sup>®</sup> II Battery Performance

### **Urgent: Medical Device Correction**

Dated July 2011

You entered the following serial number:  
NGV423092H

### **This pump is included in the population of devices containing a battery manufactured on or after March 17, 2005**

A small percentage of these pumps may exhibit Low Battery Reset, premature Elective Replacement Indicator (ERI), or premature End of Service (EOS). Additionally, for a small percentage of these pumps, the minimum timeframe of 90 days between ERI and EOS may be reduced due to this battery issue.

Refer to Medtronic's physician letter dated July 2011 for important safety information and patient management recommendations.

To obtain accurate and reliable results, you must enter the serial number accurately. If you have questions about the results of your search or the use of this web page, the following resources are available:

**Within the US:** Please contact your local Medtronic representative, or Medtronic Neuromodulation Technical Services at 1-800-707-0933

**Outside the US:** Please contact your local Medtronic representative