

Session Data Report

8840 SynchroMed II B
Programming Session

Patient Information

First Name

CARL

Last Name

ERICSON

Pump Serial Number

NGV423092H

Pump Model Number

8637-40 40 mL

Cal. Constant

120

Catheter Information

Catheter Model

8709

Volume per cm

0.0022 mL/cm

Length

89.0 cm

Removed

4.0 cm

Implanted

85.0 cm

Total Catheter Volume

0.187 mL

Current Pump Settings

Examined At

02/01/2012 MM/DD/YYYY 14:40

Last Change

11/16/2011 MM/DD/YYYY 14:11

Infusion Drug

Drug

Concentration

1 BACLOFEN

700.0 µg/mL

Infusion Mode

Simple Continuous

Drug

Dose Per Day

1 BACLOFEN

300.06 µg/day

Estimated ERI

44 Months

Reservoir Volume

7.0 mL

Low Reservoir Alarm Volume

2.0 mL

Refill Interval

11 days

Low Reservoir Alarm Date

02/12/2012 MM/DD/YYYY

Critical Alarm Interval

01:00 h:m

Non-Critical Alarm Interval

01:00 h:m

Pump Status After Update

Last Change

02/01/2012 MM/DD/YYYY 14:41

Infusion Drug

Drug

Concentration

1 BACLOFEN

700.0 µg/mL

Infusion Mode

Daily Dose Change...

0 %

Simple Continuous

Drug

Dose Per Day

1 BACLOFEN

300.06 µg/day

Estimated ERI

44 Months

Reservoir Volume

40.0 mL

Low Reservoir Alarm Volume

2.0 mL

Refill Interval

88 days

Low Reservoir Alarm Date

04/29/2012 MM/DD/YYYY

Critical Alarm Interval

01:00 h:m

Non-Critical Alarm Interval

01:00 h:m

Logs

Logs have not been read