

Session Data Report

8840 SynchroMed II B  
Programming Session

Patient Information

First Name

CARL

Last Name

ERICSON

Pump Serial Number

NGV423092H

Pump Model Number

8637-40 40 mL

Cal. Constant

120

Catheter Information

Catheter Model

8709

Volume per cm

0.0022 mL/cm

Length

89.0 cm

Removed

4.0 cm

Implanted

85.0 cm

Total Catheter Volume

0.187 mL

Current Pump Settings

Examined At

03/04/2011 MM/DD/YYYY 11:33

Last Change

12/10/2010 MM/DD/YYYY 11:01

Infusion Drug

Drug

Concentration

1 BACLOFEN

700.0 µg/mL

Infusion Mode

Simple Continuous

Drug

Dose Per Day

1 BACLOFEN

275.10 µg/day

Estimated ERI

55 Months

Reservoir Volume

7.0 mL

Low Reservoir Alarm Volume

5.0 mL

Refill Interval

5 days

Low Reservoir Alarm Date

03/09/2011 MM/DD/YYYY

**GREEN VALLEY DRUGS**

TOLL FREE (888)433-6474

1850 WHITNEY MESA DR. STE 180 HENDERSON NV 89014

DEA#: BG8108107

Ph#(702)564-2079

Fax: ( )948-6820

PIC: SCOT SILBER

Rx#: 479660

ERICSON, CARL

PHAM, TAI

4250 ARVILLE ST #373

02/25/2011

6950 W DESERT INN RD

LAS VEGAS NV 89103

LAS VEGAS NV 89117

Ph: (702)722-2658

FOR PUMP USE-----

(BACLOFEN 700MCG/ML)40ML

Qty: 40 BACLOFEN INTRATHECAL SOLN

Discard by: 05/25/2011 RPh/Tech: KL/MS Mfg:

No Refills. Authorization Required.

Pump Status After Update

Last Change

03/04/2011 MM/DD/YYYY 12:05

Infusion Drug

Drug

Concentration

1 BACLOFEN

700.0 µg/mL

Infusion Mode

Daily Dose Change

0 %

Simple Continuous

Drug

Dose Per Day

1 BACLOFEN

275.10 µg/day

Estimated ERI

55 Months

Reservoir Volume

40.0 mL

Low Reservoir Alarm Volume

3.0 mL

Refill Interval

94 days

Low Reservoir Alarm Date

06/06/2011 MM/DD/YYYY

Critical Alarm Interval

01:00 h:m

Non-Critical Alarm Interval

01:00 h:m

Logs